



WebMRE®

Save time. Capture manual test results with WebMRE.

A point of care (POC) testing system would not be complete without the ability to capture results for patient and liquid QC tests performed manually by those outside the clinical lab. WebMRE is a web-based, add-on module that fully integrates with TELCOR QML® to capture manual test results. That means all results – whether from connected devices or manually entered – are available within QML and interfaced to the LIS/EMR for consistent and consolidated point of care result reporting.

WebMRE uses the existing functions of QML. It allows the laboratory to define unlimited panel and test configurations, result options for each test, and flags for abnormal results. WebMRE uses Import and other QML functions to manage manual test results and the operators who perform them, just like it does when receiving automated test results. It also uses the existing result interface, whether EDI or script, to report these results to the LIS/EMR. Not only does this provide consolidated and consistent POC reporting, but it also provides the opportunity to bill for manual tests as a byproduct of their result entry.

WebMRE Features

Manual Result Entry

* Indicates a required field

Acute/IC/Lot Level 123456

Patient Patient

Entered By

Collection Date/Time 7/11/2022 14:37

Result Type * FIN

Acute/IC/Lot Level 123456

Collection Date * 7/11/2022

Collection Time * (Military Time - HHMM) 1437

Performed By *

Ordering Provider

Provider ID

Provider Last Name

Provider First Name

Provider Middle Initial

Test Configuration

Test Name

Lot

Device Code

Sample Type

Performed Date * (Military Time - HHMM) 1437

Comment 1

Comment 2

Comment 3

Test

Enter Pregnancy * P

Enter Pregnancy Interval QC

Verify patient ID through an ADT interface or select QC lot/level for QC results.

Enter patient and liquid QC results.

Define an unlimited number of panels, which are only accessible to operators certified to perform the panels.

Capture consistent data for each test panel.

Tests defined for each panel including three columns for the test, flag, and result values.





Everything, Everywhere Point of Care™

Save time. Improve operator management and compliance.

WebMRE not only improves your workflow, but also plays a critical role in helping you maintain compliance for your point of care program. It provides the tools you need to manage operator certifications, as well as track patient and QC results manually entered via WebMRE or electronically received from devices.

WebMRE: A winning, time-saving solution.

Laboratory
Create unlimited
customizations



- Define unlimited panels, tests within panels, and result options for each test
- Capture patient and QC results as discrete data elements
- Assign operators only to tests for which they are certified
- Access real-time result exceptions

Charting
Save time &
minimize errors



- Access WebMRE from existing PCs on your intranet
- Minimize key strokes through drop-down lists
- Give access to certified tests for only the operator signed on
- Gain result flagging for both numeric and alpha results

Improve
compliance



- Reduce errors
- Consolidate all POC results for clinician access
- Provide consistent documentation
- Capture automated electronic charges

"Before, I had no way of tracking competency with manual tests. It was a huge benefit for me to use WebMRE to track that. WebMRE helps keep everyone in line as best as we can do."

– Point of Care Coordinator,
Children's Hospitals and
Clinics of Minnesota



REACH FOR MORE >>>

Copyright 2026 TELCOR Inc. All rights reserved.

866-489-1207

telcor@telcor.com

