


# Discover why more than 2,700 facilities choose TELCOR for their POC Solution



## > UNATTENDED INTERFACES WITH UNLIMITED SUPPORT



TELCOR QML® is designed to provide a solution without limitations, supporting the management and reporting requirements of the growing point of care industry and connecting all device types—regardless of manufacturer, number of devices, or number of facilities.

QML is not only the most comprehensive point of care solution, it also connects more device types and LIS/EMR systems than anyone in the industry. This provides our customers with unparalleled point of care connectivity so you can connect “everything, everywhere, point of care®.”

## > EXTENSIVE FUNCTIONALITY & TOOLS

- ◇ Automated Alerts – Send interface alerts via text or email.
- ◇ Result Exception Management – Configured by flags, ranges, or missing required data.
- ◇ Operator Management – Full competency tracking including interfaces to LMS with automatic recertification options.
- ◇ Device Management – Bi-directional communication and full device configurations for many.
- ◇ Metrics and Statistics – Utilization, QC Levey Jennings and Linearity, TAT, and more.
- ◇ Unlimited Data Mining – Unlimited filter options and data views with print, fax, PDF, email, or export options.
- ◇ Security – Compliant with HIPAA, HITECH; access by role and area from anywhere.
- ◇ Manual Results – Unlimited configuration and charting of patient and QC through WebMRE®.

# CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
Abbott Core Laboratory	CELL-DYN Emerald	Uni-directional
Abbott Diabetes Care	Abbott Glucose Direct	Bi-directional full
Abbott Diabetes Care	Abbott Glucose PWeb	Bi-directional
Abbott Diabetes Care	Abbott Glucose UniPOC	Bi-directional
Abbott Diabetes Care	Abbott Ketone Direct	Bi-directional full
Abbott Diabetes Care	Abbott Ketone PWeb	Bi-directional
Abbott Diabetes Care	Abbott Ketone UniPOC	Bi-directional
Abbott Point of Care	Afinion	Bi-directional full
Abbott Point of Care	Alere i	Uni-directional
Abbott Point of Care	i-STAT Alinity	Bi-directional
Abbott Point of Care	i-STAT CDS5	Uni-directional
Abbott Point of Care	i-STAT DE	Bi-directional
Abbott Point of Care	i-STAT Info HQ	Bi-directional
Abbott Point of Care	i-STAT PWeb	Bi-directional
Abbott Point of Care	i-STAT UniPOC	Bi-directional
Abbott Point of Care	Piccolo xpress w Ethernet	Uni-directional
Abbott Point of Care	RALS	Bi-directional
Abbott Rapid Dx	ID NOW	Bi-directional full
Acelis Connected Health	CoagClinic	Uni-directional
Advanced Instruments	Osmo1	Uni-directional
Advanced Instruments	OsmoPRO**	Uni-directional
ALCOR Scientific	iSED/iSED ELITE**	Uni-directional
ALCOR Scientific	miniISED**	Uni-directional
Alfa Wassermann	ACE Axcel	Uni-directional
ARKRAY	Aution Eleven	Uni-directional
Beckman Coulter	Access 2 Immunoassay System**	Bi-directional
Beckman Coulter	Coulter AcT diff, 2, 8, 10, 5diff AL	Uni-directional
Beckman Coulter	Coulter DxH 500/520	Uni-directional
Becton, Dickinson and Company	Veritor Plus	Uni-directional
Binx Health	binx io	Uni-directional
bioMérieux	BIOFIRE SPOTFIRE	Bi-directional full
bioMérieux	FilmArray2.0/Torch	Uni-directional
Bio-Rad Laboratories, Inc.	in2it AIC	Uni-directional
Cepheid	GeneXpert DX	Uni-directional
Cepheid	GeneXpert Xpress	Bi-directional
Clarity Diagnostics, LLC	Clarity Platinum	Uni-directional
Clarity Diagnostics, LLC	McKesson 120	Uni-directional
Clarity Diagnostics, LLC	Urocheck 120	Uni-directional
Clinical Diagnostic Solutions, Inc.	Medonic M Series	Uni-directional
CoaguSense, Inc.	Coag-Sense	Bi-directional
CTMI Cole-Taylor Marketing, Inc.	CT-120	Uni-directional
Drager Medical, Inc.	JM-105	Uni-directional
Drucker Diagnostics	QBC STAR**	Uni-directional
EKF Life Sciences	SIRRUS Clinical Chemistry Analyzer**	Bi-directional
Haemonetics Corporation	TEG 5000	Uni-directional
Haemonetics Corporation	TEG 6s/Haemocommunicator	Bi-directional
Helena Laboratories Corp.	Actalyke XL	Uni-directional
Helena Laboratories Corp.	Cascade POC	Uni-directional
HemoCue America	201 DM Glu	Bi-directional full
HemoCue America	201 DM Hb	Bi-directional full
HemoCue America	201 DM Software	Uni-directional
HemoSonics	Quanta Hemostasis System	Bi-directional
Hemosure	Accu-Reader A100**	Uni-directional
Horiba	ABX Micros 60	Uni-directional
Horiba	ABX Pentra 60C+	Uni-directional
Horiba	ABX Pentra C400**	Uni-directional
iLine Microsystems	microINR**	Uni-directional
Intoximeters, Inc.	IntoxNet	Uni-directional

# CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
LifeHealth, LLC	IRMA LH**	Bi-directional
LifeHealth, LLC	IRMA TruPoint	Uni-directional
LSI Medience Corporation	PATHFAST	Bi-directional
LumiraDx	LumiraDx	Bi-directional full
M Dialysis AB	ISCUSflex	Uni-directional
Magellan Diagnostics	LeadCare Plus	Uni-directional
MEDTOX Scientific, Inc.	MEDTOXScan	Uni-directional
Medtronic	ACT II/InterACT	Uni-directional
Medtronic	ACT Plus	Uni-directional
Medtronic	HMS Plus	Uni-directional
Meridian Bioscience	BreathID**	Uni-directional
Nova Biomedical	NovaNet Prime/Plus	Bi-directional
Nova Biomedical	NovaNet StatSensor	Bi-directional
Nova Biomedical	NovaNet StatStrip Glucose	Bi-directional
Nova Biomedical	NovaNet StatStrip Ketone	Bi-directional
Nova Biomedical	NovaNet StatStrip Lactate	Bi-directional
Nova Biomedical	NovaNet StatStrip MTE	Bi-directional
Nova Biomedical	Stat Profile pHox Ultra	Bi-directional
Nova Biomedical	Stat Profile Prime	Bi-directional
NOWDiagnostics	ADEXUSDx	Bi-directional
OPTI Medical Systems, Inc.	OPTI CCA/CCA-R	Uni-directional
OPTI Medical Systems, Inc.	OPTI CCA-TS	Uni-directional
OPTI Medical Systems, Inc.	OPTI CCA-TS2	Uni-directional
Philips Healthcare	BiliChek	Uni-directional
PixCell Medical	HemoScreen	Bi-directional
Qorvo	Omnia	Uni-directional
Quidel	Savanna*	Bi-directional
Quidel	Sofia/Sofia 2	Bi-directional
Quidel	Solana	Uni-directional
Quidel	Triage via Census	Uni-directional
Quidel	Triage via TMSI	Bi-directional
Radiometer America	ABL700s	Bi-directional
Radiometer America	ABL77	Bi-directional
Radiometer America	ABL80	Bi-directional
Radiometer America	ABL80 FLEX	Bi-directional
Radiometer America	ABL800s	Bi-directional
Radiometer America	ABL90 FLEX/PLUS	Bi-directional
Radiometer America	AQUIRE Blood Gas	Bi-directional
Radiometer America	AQUIRE Hemocue	Uni-directional
Radiometer America	Radiance	Bi-directional
Response Biomedical	RAMP	Bi-directional
Response Biomedical	RAMP 200	Bi-directional
Roche Diagnostics	CoaguChek Pro II**	Bi-directional full
Roche Diagnostics	CoaguChek XS Plus	Bi-directional full
Roche Diagnostics	CoaguChek XS Plus w/o QCM	Bi-directional
Roche Diagnostics	CoaguChek XS Plus/Pro w/o QCM	Bi-directional
Roche Diagnostics	CoaguChek XS Pro	Bi-directional full
Roche Diagnostics	CoaguChek XS Pro w/o QCM	Bi-directional
Roche Diagnostics	cobas b 123	Uni-directional
Roche Diagnostics	cobas b 221	Uni-directional
Roche Diagnostics	cobas bge link	Bi-directional
Roche Diagnostics	Inform Glu via navify POC Operations	Bi-directional
Roche Diagnostics	Inform Glu via RALS	Uni-directional
Roche Diagnostics	Inform OTE via navify POC Operations	Bi-directional
Roche Diagnostics	Liat	Bi-directional full
Roche Diagnostics	OMNI 1-9	Uni-directional
Roche Diagnostics	Urisys 1100	Bi-directional

\*On our development schedule. \*\*Evaluation completed.

# CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
Siemens Healthcare Diagnostics	CLINITEK Connect	Bi-directional
Siemens Healthcare Diagnostics	CLINITEK Status	Uni-directional
Siemens Healthcare Diagnostics	DCA Vantage	Bi-directional
Siemens Healthcare Diagnostics	epoc direct	Bi-directional full
Siemens Healthcare Diagnostics	epoc via EDM	Bi-directional
Siemens Healthcare Diagnostics	epoc via POCcelerator SE	Bi-directional
Siemens Healthcare Diagnostics	IMMULITE 1000	Uni-directional
Siemens Healthcare Diagnostics	RAPIDComm	Bi-directional
Siemens Healthcare Diagnostics	RAPIDComm BG	Bi-directional
Siemens Healthcare Diagnostics	RAPIDComm DCA	Bi-directional
Siemens Healthcare Diagnostics	RAPIDComm UA	Bi-directional
Siemens Healthcare Diagnostics	RAPIDPoint	Bi-directional
Siemens Healthcare Diagnostics	Stratus CS	Uni-directional
Siemens Healthcare Diagnostics	UniPOC	Bi-directional
Sight Diagnostics LTD	OLO	Uni-directional
Stago	STA Compact Max	Uni-directional
Sysmex America, Inc.	poch-100i	Uni-directional
Sysmex America, Inc.	Sysmex XN/XN-L	Uni-directional
Sysmex America, Inc.	Sysmex XP	Uni-directional
Sysmex America, Inc.	Sysmex XS	Uni-directional
Sysmex America, Inc.	Sysmex XS 1000i	Uni-directional
Sysmex America, Inc.	Sysmex XT	Uni-directional
Universal Biosensors	Xprecia Prime	Bi-directional
Werfen North America	AVOX 1000E	Uni-directional
Werfen North America	AVOX 4000	Uni-directional
Werfen North America	GEM 3000	Bi-directional
Werfen North America	GEM 3500	Bi-directional
Werfen North America	GEM 4000	Bi-directional
Werfen North America	GEM 5000	Bi-directional
Werfen North America	GEM 7000	Bi-directional
Werfen North America	GEM OPL	Uni-directional
Werfen North America	GEMWeb Plus	Bi-directional
Werfen North America	GEMWebPlus ChemSTAT	Bi-directional
Werfen North America	GH100	Bi-directional full
Werfen North America	ROTEM Delta	Bi-directional
Werfen North America	ROTEM Delta via GWP	Bi-directional
Werfen North America	ROTEM Sigma via GWP	Bi-directional
Werfen North America	Sig Elite	Bi-directional
Werfen North America	VerifyNow	Uni-directional

**\*On our development schedule.**

**\*\*Evaluation completed**

**Uni-directional:** may include patient, AQC, and EQC results to QML.

**Bi-directional:** uni-directional and may include date/time, some device configurations, operators, patient lists, reagents, device events.

**Bi-directional full:** uni-directional, bi-directional, and may include full device configurations, firmware/software updates, languages.

Functionality is dependent on the device interface specifications provided by the manufacturer.

All trademarks, service marks, trade names, trade dress, product names and logos appearing on this document are the property of their respective owners.

