

Discover why more than 2,600 facilities choose TELCOR for their POC Solution



▶ UNATTENDED INTERFACES WITH UNLIMITED SUPPORT



TELCOR QML® is designed to provide a solution without limitations, supporting the management and reporting requirements of the growing point of care industry and connecting all device types—regardless of manufacturer, number of devices, or number of facilities.

QML is not only the most comprehensive point of care solution, it also connects more device types and LIS/EMR systems than anyone in the industry. This provides our customers with unparalleled point of care connectivity so you can connect “everything, everywhere, point of care™.”

▶ EXTENSIVE FUNCTIONALITY & TOOLS

- ◊ Automated Alerts – Send interface alerts via text or email
- ◊ Result Exception Management – Configured by flags, ranges, or missing required data
- ◊ Operator Management – Full competency tracking including interfaces to LMS with automatic recertification options
- ◊ Device Management – Bi-directional communication and full device configurations for many
- ◊ Metrics and Statistics – Utilization, QC Levey Jennings and Linearity, TAT, and more
- ◊ Unlimited Data Mining – Unlimited filter options and data views with print, fax, PDF, email, or export options
- ◊ Security – Compliant with HIPAA, HITECH; access by role and area from anywhere
- ◊ Manual Results – Unlimited configuration and charting of patient and QC through WebMRE®

▶ *Discover more at TELCOR.com or call 866-489-1207.*

CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
Abbott Core Laboratory	Cell-DynEmerald	Uni-directional
Abbott Diabetes Care	Abbott Glucose Direct	Full bi-directional
Abbott Diabetes Care	Abbott Glucose PWeb	Limited bi-directional
Abbott Diabetes Care	Abbott Glucose UniPOC	Limited bi-directional
Abbott Diabetes Care	Abbott Ketone Direct	Full bi-directional
Abbott Diabetes Care	Abbott Ketone PWeb	Limited bi-directional
Abbott Diabetes Care	Abbott Ketone UniPOC	Limited bi-directional
Abbott Point of Care	Afinion	Full bi-directional
Abbott Point of Care	Alere i	Uni-directional
Abbott Point of Care	i-STAT Alinity	Limited bi-directional
Abbott Point of Care	i-STAT CDS5	Uni-directional
Abbott Point of Care	i-STAT DE	Limited bi-directional
Abbott Point of Care	i-STAT Info HQ	Limited bi-directional
Abbott Point of Care	i-STAT PWeb	Limited bi-directional
Abbott Point of Care	i-STAT UniPOC	Limited bi-directional
Abbott Point of Care	Piccolo xpress w Ethernet	Uni-directional
Abbott Point of Care	RALS	Limited bi-directional
Abbott Rapid Dx	ID NOW	Full bi-directional
Acelis Connected Health	CoagClinic	Uni-directional
Advanced Instruments	Osmo1	Uni-directional
Advanced Instruments	OsmoPRO**	Uni-directional
Alfa Wassermann	ACE Axcel	Uni-directional
ARKRAY	Aution Eleven	Uni-directional
Beckman Coulter	Access 2 Immunoassay System*	Limited bi-directional
Beckman Coulter	Coulter AcT diff, 2, 8, 10, 5diff AL	Uni-directional
Beckman Coulter	Coulter DxH 500/520	Uni-directional
Becton, Dickinson and Company	Veritor Plus	Uni-directional
Binx Health	binx io*	Uni-directional
bioMérieux	BIOFIRE SPOTFIRE	Full bi-directional
Bio-Rad Laboratories, Inc.	in2it AIC	Uni-directional
Cepheid	GeneXpert DX	Uni-directional
Cepheid	GeneXpert Xpress	Limited bi-directional
Clarity Diagnostics, LLC	Clarity Platinum	Uni-directional
Clarity Diagnostics, LLC	McKesson 120	Uni-directional
Clarity Diagnostics, LLC	Urocheck 120	Uni-directional
Clinical Diagnostic Solutions, Inc.	Medonic M Series	Uni-directional
CoaguSense, Inc.	Coag-Sense	Limited bi-directional
CTMI Cole-Taylor Marketing, Inc.	CT-120	Uni-directional
Drager Medical, Inc.	JM-105	Uni-directional
Drucker Diagnostics	QBC STAR**	Uni-directional
EKF Life Sciences	SIRRUS Clinical Chemistry Analyzer**	Limited bi-directional
Haemonetics Corporation	TEG 5000	Uni-directional
Haemonetics Corporation	TEG 6s/Haemocommunicator	Limited bi-directional
Helena Laboratories Corp.	Actalyke XL	Uni-directional
Helena Laboratories Corp.	Cascade POC	Uni-directional
HemoCue America	201 DM Glu	Full bi-directional
HemoCue America	201 DM Hb	Full bi-directional
HemoCue America	201 DM Software	Uni-directional
HemoSonics	Quantra Hemostasis System	Limited bi-directional
Hemosure	Accu-Reader A100**	Uni-directional
Horiba	ABX Micros 60	Uni-directional
Horiba	ABX Pentra 60C+	Uni-directional
Horiba	ABX Pentra C400**	Uni-directional
iLine Microsystems	microINR**	Uni-directional
Intoximeters, Inc.	IntoxNet	Uni-directional
LifeHealth, LLC	IRMA LH**	Limited bi-directional
LifeHealth, LLC	IRMA TruPoint	Uni-directional
LifeScan, a Johnson & Johnson	SureStep Pro/Flex	Limited bi-directional
LSI Medience Corporation	PATHFAST	Limited bi-directional

CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
LumiraDx	LumiraDx	Full bi-directional
M Dialysis AB	ISCUflex	Uni-directional
Magellan Diagnostics	LeadCare Plus	Uni-directional
MEDTOX Scientific, Inc.	MEDTOXScan	Uni-directional
Medtronic	ACT II/InterACT	Uni-directional
Medtronic	ACT Plus	Uni-directional
Medtronic	HMS Plus	Uni-directional
Meridian bioscience	BreathID**	Uni-directional
Nova Biomedical	NovaNet Prime/Plus	Limited bi-directional
Nova Biomedical	NovaNet StatSensor	Limited bi-directional
Nova Biomedical	NovaNet StatStrip Glucose	Limited bi-directional
Nova Biomedical	NovaNet StatStrip Ketone	Limited bi-directional
Nova Biomedical	NovaNet StatStrip Lactate	Limited bi-directional
Nova Biomedical	NovaNet StatStrip MTE	Limited bi-directional
Nova Biomedical	PDM	Limited bi-directional
Nova Biomedical	Stat Profile pHox Ultra	Limited bi-directional
Nova Biomedical	Stat Profile Prime	Limited bi-directional
NOWDiagnostics	ADEXUSdx	Limited bi-directional
OPTI Medical Systems, Inc.	OPTI CCA/CCA-R	Uni-directional
OPTI Medical Systems, Inc.	OPTI CCA-TS	Uni-directional
OPTI Medical Systems, Inc.	OPTI CCA-TS2	Uni-directional
Philips Healthcare	BiliChek	Uni-directional
PixCell Medical	HemoScreen	Limited bi-directional
Qorvo	Omnia	Uni-directional
Quidel	Savanna*	Limited bi-directional
Quidel	Sofia/Sofia 2	Limited bi-directional
Quidel	Triage via Census	Uni-directional
Quidel	Triage via TMSI	Limited bi-directional
Radiometer America	ABL700s	Limited bi-directional
Radiometer America	ABL77	Limited bi-directional
Radiometer America	ABL80	Limited bi-directional
Radiometer America	ABL80 FLEX	Limited bi-directional
Radiometer America	ABL800s	Limited bi-directional
Radiometer America	ABL90 FLEX/PLUS	Limited bi-directional
Radiometer America	AQT90 FLEX	Limited bi-directional
Radiometer America	AQUIRE Blood Gas	Limited bi-directional
Radiometer America	AQUIRE Hemocue	Uni-directional
Radiometer America	Radiance	Limited bi-directional
Response Biomedical	RAMP	Limited bi-directional
Response Biomedical	RAMP 200	Limited bi-directional
Roche Diagnostics	CoaguChek Pro II**	Full bi-directional
Roche Diagnostics	CoaguChek XS Plus	Full bi-directional
Roche Diagnostics	CoaguChek XS Plus limited	Limited bi-directional
Roche Diagnostics	CoaguChek XS Plus/Pro limited	Limited bi-directional
Roche Diagnostics	CoaguChek XS Pro	Full bi-directional
Roche Diagnostics	CoaguChek XS Pro limited	Limited bi-directional
Roche Diagnostics	cobas b 123	Uni-directional
Roche Diagnostics	cobas b 221	Uni-directional
Roche Diagnostics	cobas bge link	Limited bi-directional
Roche Diagnostics	Inform Glu via infinity POC	Limited bi-directional
Roche Diagnostics	Inform Glu via RALS	Uni-directional
Roche Diagnostics	Inform OTE via infinity POC	Limited bi-directional
Roche Diagnostics	Liat	Full bi-directional
Roche Diagnostics	OMNI 1-9	Uni-directional
Roche Diagnostics	Urisys 1100	Limited bi-directional
Siemens Healthcare Diagnostics	Atellica DCA*	Limited bi-directional
Siemens Healthcare Diagnostics	CLINITEK Connect	Limited bi-directional

*On our development schedule. **Evaluation completed.

CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
Siemens Healthcare Diagnostics	CLINITEK Status	Uni-directional
Siemens Healthcare Diagnostics	DCA Vantage	Limited bi-directional
Siemens Healthcare Diagnostics	epoc direct	Full bi-directional
Siemens Healthcare Diagnostics	epoc via EDM	Limited bi-directional
Siemens Healthcare Diagnostics	epoc via POCcelerator SE	Limited bi-directional
Siemens Healthcare Diagnostics	IMMULITE 1000	Uni-directional
Siemens Healthcare Diagnostics	RAPIDComm	Limited bi-directional
Siemens Healthcare Diagnostics	RAPIDComm BG	Limited bi-directional
Siemens Healthcare Diagnostics	RAPIDComm DCA	Limited bi-directional
Siemens Healthcare Diagnostics	RAPIDComm UA	Limited bi-directional
Siemens Healthcare Diagnostics	RAPIDPoint	Limited bi-directional
Siemens Healthcare Diagnostics	Stratus CS	Uni-directional
Siemens Healthcare Diagnostics	UniPOC	Limited bi-directional
Sight Diagnostics LTD	OLO	Uni-directional
Stago	STA Compact Max	Uni-directional
Sysmex America, Inc.	pocH-100i	Uni-directional
Sysmex America, Inc.	Sysmex XN/XN-L	Uni-directional
Sysmex America, Inc.	Sysmex XP	Uni-directional
Sysmex America, Inc.	Sysmex XS	Uni-directional
Sysmex America, Inc.	Sysmex XS 1000i	Uni-directional
Sysmex America, Inc.	Sysmex XT	Uni-directional
Werfen North America	AVOX 1000E	Uni-directional
Werfen North America	AVOX 4000	Uni-directional
Werfen North America	GEM 3000	Limited bi-directional
Werfen North America	GEM 3500	Limited bi-directional
Werfen North America	GEM 4000	Limited bi-directional
Werfen North America	GEM 5000	Limited bi-directional
Werfen North America	GEM OPL	Uni-directional
Werfen North America	GEMWeb Plus	Limited bi-directional
Werfen North America	GEMWebPlus ChemSTAT	Limited bi-directional
Werfen North America	GH100	Full bi-directional
Werfen North America	ROTEM Delta	Limited bi-directional
Werfen North America	ROTEM Delta via GWP	Limited bi-directional
Werfen North America	ROTEM Sigma via GWP	Limited bi-directional
Werfen North America	Sig Elite	Limited bi-directional
Werfen North America	VerifyNow	Uni-directional

***On our development schedule.**

****Evaluation completed**

Uni-directional: may include patient, AQC, and EQC results to QML

Limited bi-directional: uni-directional plus may include date/time, operators, patient lists, reagents.

Full bi-directional: uni-directional, limited bi-directional, and may include device configurations, firmware/software updates, languages.

All trademarks, service marks, trade names, trade dress, product names and logos appearing on this document are the property of their respective owners.

About Us

TELCOR provides connectivity and data management to more POC device types and LIS/EMR systems than any provider in the industry. We are the only POC middleware vendor with no diagnostic or LIS company ownership allowing us to connect any manufacturer's devices. With more than 10,000 device types at more than 2,500 hospitals and thousands of ambulatory sites connected to TELCOR QML, facilities are able to collect results electronically thus saving time and improving data quality. By focusing on software, we are not distracted with selling devices or reagents. We positively provide the best solutions, support, and service in the industry.

Get in Touch



7101 A Street
Lincoln, Nebraska 68510



sales@telcor.com | telcor@telcor.com



866-489-1207
402-489-1207