

Discover why more than 2,500 facilities choose TELCOR for their POC Solution



▶ UNATTENDED INTERFACES WITH UNLIMITED SUPPORT



TELCOR QML® is designed to provide a solution without limitations, supporting the management and reporting requirements of the growing point of care industry and connecting all device types—regardless of manufacturer, number of devices, or number of facilities.

QML is not only the most comprehensive point of care solution, it also connects more device types and LIS/EMR systems than anyone in the industry. This provides our customers with unparalleled point of care connectivity so you can connect “everything, everywhere, point of care™.”

▶ EXTENSIVE FUNCTIONALITY & TOOLS

- ◇ Automated Alerts – Send interface alerts via text or email
- ◇ Result Exception Management – Configured by flags, ranges, or missing required data
- ◇ Operator Management – Full competency tracking including interfaces to LMS with automatic recertification options
- ◇ Device Management – Bi-directional communication and full device configurations for many
- ◇ Metrics and Statistics – Utilization, QC Levey Jennings and Linearity, TAT, and more
- ◇ Unlimited Data Mining – Unlimited filter options and data views with print, fax, PDF, email, or export options
- ◇ Security – Compliant with HIPAA, HITECH; access by role and area from anywhere
- ◇ Manual Results – Unlimited configuration and charting of patient and QC through WebMRE®

▶ *Discover more at TELCOR.com or call 866-489-1207.*

CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
Abbott Core Laboratory	Cell-DynEmerald	Uni-directional
Abbott Diabetes Care	Abbott Glucose Direct	Full bi-directional
Abbott Diabetes Care	Abbott Glucose PWeb, UniPOC	Limited bi-directional
Abbott Diabetes Care	Abbott Ketone Direct	Full bi-directional
Abbott Diabetes Care	Abbott Ketone PWeb, UniPOC	Limited bi-directional
Abbott Point of Care	Afinion	Full bi-directional
Abbott Point of Care	Alere i	Uni-directional
Abbott Point of Care	ID NOW	Full bi-directional
Abbott Point of Care	i-STAT Alinity	Limited bi-directional
Abbott Point of Care	i-STAT CDS5	Uni-directional
Abbott Point of Care	i-STAT DE	Limited bi-directional
Abbott Point of Care	i-STAT Info HQ	Uni-directional
Abbott Point of Care	i-STAT PWeb, UniPOC	Limited bi-directional
Abbott Point of Care	Piccolo xpress	Uni-directional
Abbott Point of Care	RALS	Limited bi-directional
Acelis Connected Health	CoagClinic	Uni-directional
Advanced Instruments	Osmol	Uni-directional
Alfa Wassermann	ACE Axcel	Uni-directional
ARKRAY	Aution Eleven	Uni-directional
Beckman Coulter	Coulter AcT diff, 2, 8, 10, 5diff AL	Uni-directional
Beckman Coulter	Coulter DxH 500/520	Uni-directional
Becton, Dickinson and Company	Veritor Plus	Uni-directional
Bio-Rad Laboratories, Inc.	in2it AIC	Uni-directional
Cepheid	GeneXpert DX	Uni-directional
Cepheid	GeneXpert Xpress	Limited bi-directional
Clarity Diagnostics, LLC	McKesson 120	Uni-directional
Clarity Diagnostics, LLC	Urocheck 120	Uni-directional
Clinical Diagnostic Solutions, Inc.	Medonic M Series	Uni-directional
CoaguSense, Inc.	Coag-Sense	Limited bi-directional
CTMI Cole-Taylor Marketing, Inc.	CT-120	Uni-directional
Drager Medical, Inc.	JM-105	Uni-directional
EKF Life Sciences	SIRRUS Clinical Chemistry Analyzer	Limited bi-directional
Haemonetics Corporation	TEG 5000	Uni-directional
Haemonetics Corporation	TEG 6s/Haemocommunicator	Limited bi-directional
Helena Laboratories Corp.	Actalyke XL	Uni-directional
Helena Laboratories Corp.	Cascade POC	Uni-directional
HemoCue America	HemoCue 201 DM	Uni-directional
HemoCue America	HemoCue Direct*	Full bi-directional
HemoSonics	Quantra QPlus System	Limited bi-directional
Horiba	ABX Micros 60, 60C+	Uni-directional
Intoximeters, Inc.	IntoxNet	Uni-directional
LifeHealth, LLC	IRMA TruPoint	Uni-directional
LifeScan	SureStep Pro/Flex	Limited bi-directional
LSI Medience Corporation	PATHFAST	Limited bi-directional
LumiraDx	LumiraDx*	Full bi-directional
M Dialysis AB	ISCUSflex	Uni-directional
Magellan Diagnostics	LeadCare Plus	Uni-directional
MEDTOX Scientific, Inc.	MEDTOXScan	Uni-directional
Medtronic	ACT II/InterACT/ACT Plus	Uni-directional
Medtronic	HMS Plus	Uni-directional
Nova Biomedical	NovaNet Prime/Plus	Limited bi-directional
Nova Biomedical	NovaNet StatSensor, Ketone, Lactate	Limited bi-directional
Nova Biomedical	NovaNet StatStrip Glucose, MTE	Limited bi-directional
Nova Biomedical	PDM	Uni-directional
Nova Biomedical	Stat Profile pHox Ultra	Uni-directional
Nova Biomedical	Stat Profile Prime	Uni-directional
NOWDiagnostics	ADEXUSDx	Limited bi-directional
OPTI Medical Systems, Inc.	OPTI CCA/CCA-R	Uni-directional
OPTI Medical Systems, Inc.	OPTI CCA-TS	Uni-directional

CONNECTED DEVICE TYPES

Vendor	Device Name	Functionality
OPTI Medical Systems, Inc.	OPTI CCA-TS2	Uni-directional
Philips Healthcare	BiliChek	Uni-directional
PixCell Medical	HemoScreen*	Full bi-directional
Quidel	Sofia/Sofia 2	Limited bi-directional
Quidel	Triage via Census	Uni-directional
Quidel	Triage via TMSI	Limited bi-directional
Radiometer America	ABL700s, ABL77	Limited bi-directional
Radiometer America	ABL80, ABL80 FLEX, ABL800s	Limited bi-directional
Radiometer America	ABL90 FLEX/PLUS	Limited bi-directional
Radiometer America	AQT90 FLEX	Limited bi-directional
Radiometer America	AQURE Blood Gas	Limited bi-directional
Radiometer America	AQURE Hemocue	Uni-directional
Radiometer America	Radiance	Limited bi-directional
Response Biomedical	RAMP, RAMP 200	Limited bi-directional
Roche Diagnostics	CoaguChek XS Plus*	Full bi-directional
Roche Diagnostics	CoaguChek XS Plus/Pro	Limited bi-directional
Roche Diagnostics	CoaguChek XS Pro*	Full bi-directional
Roche Diagnostics	cobas b 123, cobas b 221	Uni-directional
Roche Diagnostics	cobas bge link	Limited bi-directional
Roche Diagnostics	Inform Glu via cobas infinity	Limited bi-directional
Roche Diagnostics	Inform OTE via cobas infinity	Limited bi-directional
Roche Diagnostics	Liat	Full bi-directional
Roche Diagnostics	OMNI 1-9	Uni-directional
Roche Diagnostics	Roche Inform via RALS	Limited bi-directional
Roche Diagnostics	Urisys 1100	Limited bi-directional
Siemens Healthcare Diagnostics	CLINITEK Connect	Limited bi-directional
Siemens Healthcare Diagnostics	CLINITEK Status	Uni-directional
Siemens Healthcare Diagnostics	DCA Vantage	Limited bi-directional
Siemens Healthcare Diagnostics	epoc direct	Full bi-directional
Siemens Healthcare Diagnostics	epoc via EDM, POCcelerator SE	Limited bi-directional
Siemens Healthcare Diagnostics	IMMULITE 1000	Uni-directional
Siemens Healthcare Diagnostics	RAPIDComm, BG, DCA, UA	Limited bi-directional
Siemens Healthcare Diagnostics	RAPIDPoint	Uni-directional
Siemens Healthcare Diagnostics	Stratus CS	Uni-directional
Siemens Healthcare Diagnostics	UniPOC	Limited bi-directional
Sight Diagnostics LTD	OLO	Uni-directional
Stago	STA Compact Max	Uni-directional
Sysmex America, Inc.	poch-100i	Uni-directional
Sysmex America, Inc.	Sysmex XN/XN-L	Uni-directional
Sysmex America, Inc.	Sysmex XP	Uni-directional
Sysmex America, Inc.	Sysmex XS, XS 1000i	Uni-directional
Sysmex America, Inc.	Sysmex XT	Uni-directional
Werfen North America	AVOX 1000E, 4000	Uni-directional
Werfen North America	GEM 3000, 3500	Uni-directional
Werfen North America	GEM 4000, 5000	Limited bi-directional
Werfen North America	GEM Hemochron 100*	Full bi-directional
Werfen North America	GEM OPL	Uni-directional
Werfen North America	GEMWeb Plus	Limited bi-directional
Werfen North America	GEMWeb Plus ChemSTAT	Limited bi-directional
Werfen North America	ROTEM	Uni-directional
Werfen North America	ROTEM via GEMWeb Plus	Limited bi-directional
Werfen North America	Sig Elite	Limited bi-directional
Werfen North America	VerifyNow	Uni-directional

***On our development schedule.**

Uni-directional: may include patient, ACQ, and EQC results to QML.

Limited bi-directional: uni-directional plus may include date/time, operators, patient lists, reagents.

Full bi-directional: uni-directional, limited bi-directional, and may include device configurations, firmware/software updates, languages.

All trademarks, service marks, trade names, trade dress, product names and logos appearing on this document are the property of their respective owners.

About Us

Get in Touch

TELCOR provides connectivity and data management to more POC device types and LIS/EMR systems than any provider in the industry. We are the only POC middleware vendor with no diagnostic or LIS company ownership allowing us to connect any manufacturer's devices. With more than 10,000 device types at more than 2,500 hospitals and thousands of ambulatory sites connected to TELCOR QML, facilities are able to collect results electronically thus saving time and improving data quality. By focusing on software, we are not distracted with selling devices or reagents. We positively provide the best solutions, support, and service in the industry.



7101 A Street
Lincoln, Nebraska 68510



sales@telcor.com | telcor@telcor.com



866-489-1207
402-489-1207